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## Instructions for use CAVI-T® Balloon

### Table of content

1.	<b>Indication</b> .....	7
2.	<b>Contra-indications</b> .....	7
3.	<b>Description</b> .....	8
3.1.	Pictograms .....	8
3.2.	Composition of the device .....	8
4.	<b>Product use</b> .....	8
4.1.	<b>Step 1:</b> Device preparation.....	8
4.2.	<b>Step 2:</b> Insertion and positioning .....	9
4.3.	<b>Step 3:</b> Inflation.....	9
4.4.	<b>Step 4:</b> Patient follow-up procedure.....	11
4.5.	Nose rinsing protocol .....	11
4.6.	<b>Step 5:</b> Deflation and removal .....	11
4.7.	Potential default situations .....	12
4.8.	<b>Step 6:</b> Product scrapping .....	12

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### 1. Indication

The CAVI-T® balloon is indicated for intranasal bleeding by compression (packing) of intranasal structures.

CAVI-T® is a class IIa medical device. Please read carefully instructions for use before procedure initiation.

### 2. Contraindications

CAVI-T® use is contraindicated in the following situations:

- Severe cranio-facial trauma and / or polytrauma
- Suspicion of benign or malignant tumor of the nasal fossa (e.g. nasopharyngeal fibroma)
- Non-adult population < 18 years

**Note:** In case of significant pain, posterior or anterior migration of the balloon, a septal deviation or anatomical variation should be considered and the patient managed with particular caution.

### 3. Description

#### 3.1. Pictograms

- Be careful
- Address of legal Manufacturer
- Product reference
- Batch number
- Read instructions for use (electronic multilingual IFU)
- Manufacturing date
- Expiry date
- Product sterilized with ethylene oxide  
Single sterile barrier system
- Do not resterilize
- Do not reuse
- Medical device that complies with requirements of 93/42/EEC Directive  
Do not use if packaging is damaged
- Store away from direct sunlight
- Latex free
- Phtalates free
- Medical device
- Unique device identification

#### 3.2. Composition of the device

CAVI-T® is divided into 3 main parts :

- A deployment guide
- An asymmetrical balloon in polyurethane (PU)
- A bidirectional valve

Usage of CAVI-T® has been validated with a Becton Dickinson 30 ml Luer Lock syringe.

### 4. Product Use

#### 4.1. **Step 1:** Device preparation

- Remove the device from its sterile packaging while respecting the indicated opening direction.
- Grab the device by the deployment guide (white part).
- Verify markers alignment on the guide.

Markers must be aligned and presented superiorly prior to introduction in order to ensure the introduction in the correct direction. Introducing the balloon in the incorrect direction may harm the patient and/or make the balloon ineffective.


If the product is used in the operating room, ensure the environment is sterile before using the device in order to reduce the risk of infection.

#### 4.2. **Step 2:** Insertion and positioning

The CAVI-T® balloon must be inserted in the nasal cavity with the dedicated deployment guide.

The device, placed horizontally (parallel to the bony palate) is pushed in the nasal cavity until the first marker is reached (black marker located in the same axis as the deployment guide).

In order to properly position the CAVI-T® balloon in the nasal fossa, it can be pushed till the second (marked with a « 1 », i.e. 1 cm from the proximal part of the balloon) or third marker (marked « 2 » and « MAX », i.e. 2 cm from the proximal part of the balloon).

Introduction of the balloon must be managed by a healthcare professional that has been trained to use the product. An incorrect positioning of the balloon may lead in rare instances to tissue squeezing, rupture of particular structures which could have serious consequences for the patient.

Do not push the balloon beyond the recommended upper limit (number « 2 » marked in black supplemented with the mark « MAX », both located on the deployment guide).

In case the CAVI-T® balloon placement is too posterior, it may be more painful for the patient, trigger an emetic reflex or lead to posterior migration (posterior to the sphenopalatin artery, towards the pharynx) upon inflation, thus compromising product effectiveness.

Indicators allow to position the CAVI-T® balloon adequately in the nasal cavity prior to inflation.

**Note:** optionally, the entire CAVI-T® device can be dipped into saline solution or sterile water in order to improve product gliding and improve patient comfort upon insertion and inflation. Nevertheless, this step is not mandatory if a topical anesthetic spray (e.g. xylocaine) is given to the patient prior to balloon introduction.


Respect markers printed on the device.

#### 4.3. **Step 3:** Inflation

Once CAVI-T® is correctly positioned in the nasal cavity, the balloon is inflated with air in a sequential manner with a standard Luer Lock syringe of 30 ml minimum, as described hereunder:

- Pull on the Luer Lock syringe piston until reaching 25 ml.
- Connect the syringe to the balloon valve.
- Gradually push the piston until reaching 10 ml, while controlling the patient's pain, maintaining the guide in place and exerting a slight traction on the distal portion of the tube connected to the valve to avoid posterior displacement.

This 10 ml<sup>(1)</sup> pre-inflation ensures progressive unfolding of the device.



**Note:** Pre-inflation must be performed in a controlled way. If pain occurs, stop inflation. Then, resume inflation progressively. This stepwise inflation allows the patient to get accustomed to the compression feeling.

- Once pre-inflation is finalized, remove the deployment guide and unscrew the syringe (syringe may also be left in place if maneuverability is not impacted). The valve then locks automatically.
- If needed, adjust the proximal position of the CAVI-T<sup>®</sup> balloon with the green atraumatic rounder tip located at the distal end of the rigid portion of the deployment guide while grasping the tube connected to the valve to avoid any posterior displacement.

**!** Do not use surgical forceps to reposition CAVI-T<sup>®</sup>, as this could damage the device.

In order to resume inflation, the syringe is screwed back on and the remaining air pushed until bleeding is stopped, without going above 25 ml and while maintaining the tube connected to the valve. The quantity of air injected in the balloon is constantly controlled thanks to the graduations printed on the syringe. Disconnect the syringe and verify the balloon stability by exerting a series of tractions and pressures with a (non-sharp) forceps.

<sup>1</sup>This value can change slightly according to the patient's anatomy and to the acceptable pain threshold.

**Note:** A standard Luer Lock syringe (or any other inflation device with a Luer Lock connector) is used for inflation and deflation of the device. The syringe is not included in the packaging.

**!** Carefully fix and secure the distal part of the tube on the patient's cheek with plaster (as close as possible to the nostril). Particular attention has to be paid to secure the valve in patients who underwent sinus surgery. Moreover, a more important inflation of the balloon without exceeding 25 ml, could be necessary for these patients. In some patients, overinflation can lead to a potential posterior migration of the balloon. A slight deflation often allows to correct this problem while maintaining a good level of efficacy.

#### 4.4. **Step 4:** Patient follow-up procedure

The CAVI-T<sup>®</sup> device can be used to stop intranasal bleeding in awake or non-awake patients. In both cases, a vigilance procedure is recommended for healthcare professionals in order to prevent complications, pain and discomfort, and ensure that the device works adequately.

#### Within 2 hours following balloon placement:

##### *The patient is awake:*

- 1) Perform a visual control: check product positioning and potential bleeding recurrence.
- 2) Proceed with patient questioning: pain, comfort, feeling of balloon displacement.
- 3) Verify whether the patient spits blood, which would indicate posterior bleeding recurrence.

##### *The patient is not awake:*

Only a visual control can be performed in order to check the correct positioning of the CAVI-T<sup>®</sup> balloon.

CAVI-T<sup>®</sup> is designed to remain in position in the nasal fossa for a **maximum of 72h**.


#### 4.5. Nose rinsing protocol

In order to reduce the risk of intra-balloon breathing tube obstruction, it is recommended to rinse regularly with saline or sterile water with a 10 ml syringe minimum.

Head inclined forward, the patient places the tip of the syringe (previously filled with saline or sterile water), inside the breathing tube and vigorously flushes its content. This procedure can be repeated as many times as necessary.

#### 4.6. **Step 5:** Deflation and removal

The bidirectional valve located on the CAVI-T<sup>®</sup> balloon can be used to inflate or deflate it according to needs.

CAVI-T<sup>®</sup> can be deflated to relieve the patient in case of over-inflation or to prepare for device removal.

In both cases, a standard Luer Lock syringe is connected to the valve and the balloon is emptied by pulling on the piston handle. Deflation is progressive (with a permanent control of patient pain) until the product is completely emptied of the air that it contained.

Under-inflation of the balloon may cause its posterior migration. This migration could lead to airway obstruction.

**!** Once the CAVI-T<sup>®</sup> balloon is completely emptied, it can be removed by gently pulling on the valve tubing or on the breathing tube with a (non-sharp) forceps.

#### 4.7. Potential default situations

Situation	Recommended action
The balloon cannot be inflated after the syringe has been connected.	Unscrew the syringe and firmly screw it back.
The balloon does not properly deploy once in place.	Deflate and gently remove the balloon. Perform inflation outside of the patient's nasal cavity. If control is positive, deflate the balloon, insert it in the patient's nose again and inflate it. If control is negative (e.g. leak), discard the device and pick a new one.
The balloon migrates posteriorly or is not adequately positioned in the nasal fossa.	Deflate the balloon, reposition it and re-inflate it.
The balloon cannot be removed properly from the nasal cavity (e.g. valve dysfunction).	Cut off the tubing located between the valve and the balloon for a progressive deflation. If necessary, perform a gentle traction of the balloon to bring it as close as possible to the nostril and puncture it with a sharp instrument to facilitate full extraction. Replace with a new CAVI-T <sup>®</sup> device.

#### 4.8. **Step 6:** Product scrapping

**!** The balloon should not be re-used or re-sterilized. Device reuse exposes the patient to infectious risk. Performances and sterility of the device are not insured in case of re-sterilization.

Once used, the CAVI-T<sup>®</sup> device must be discarded according to local regulation for medical devices.